# MANUAL DE USUARIO DE CONVERSOR

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Button(s)	Description
<b>Wheel</b>	Zooms camera
<b>Left</b>	Rotates camera
<b>Left double click</b>	Aligns camera
<b>Right</b>	Translates camera
<b>Right double click</b>	Centers scene
<b>Middle</b>	Zooms camera
<b>Middle double click</b>	Shows entire scene
<b>Left &amp; Middle</b>	Rotates on screen camera
<b>Shift+Left</b>	Selects
<b>Shift+Middle</b>	Zooms on region for camera
<b>Ctrl+Wheel</b>	Zooms manipulated frame
<b>Ctrl+Left</b>	Rotates manipulated frame
<b>Ctrl+Left double click</b>	Aligns frame
<b>Ctrl+Right</b>	Translates manipulated frame
<b>Ctrl+Right double click</b>	Centers frame
<b>Ctrl+Middle</b>	Zooms manipulated frame
<b>Ctrl+Left &amp; Middle</b>	Rotates on screen manipulated frame
<b>Right double click with Left pressed</b>	Resets revolve around point
<b>Left double click with Right pressed</b>	Sets revolve around point
<b>Left double click with Middle pressed</b>	Zooms on pixel
<b>Right double click with Middle pressed</b>	Zooms to fit scene

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Key(s)	Description
<b>Space</b>	Changes camera mode (revolve or fly)
<b>A</b>	Toggles the display of the world axis
<b>C</b>	Toggles camera paths display
<b>F</b>	Toggles the display of the FPS
<b>G</b>	Toggles the display of the XY grid
<b>H</b>	Opens this help window
<b>S</b>	Toggles stereo display
<b>Esc</b>	Exits program
<b>Return</b>	Starts/stops the animation
<b>Left</b>	Moves camera left
<b>Up</b>	Moves camera up
<b>Right</b>	Moves camera right
<b>Down</b>	Moves camera down
<b>Shift+?</b>	Toggles the display of the text
<b>Ctrl+S</b>	Saves a screenshot
<b>Alt+Return</b>	Toggles full screen display
Camera paths are controlled using F1..F12 (noted <i>Fx</i> below):	
<b><i>Fx</i></b>	Plays path (or resets saved position)
<b>Alt+<i>Fx</i></b>	Adds a key frame (or defines a position)
<b>Alt+<i>Fx</i>+<i>Fx</i></b>	Deletes path (or saved position)

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